

premium products - perfectly packaged



Highest quality bedding materials made from European soft woods, manufactured to ISO 9001:2008 standards exclusively for Datesand. With full batch traceability and certificates of analysis, all grades conform to GLP standards.

- Free flowing, highest quality wood granules
- Heat treated at 100°C
- Moisture content < 16%
- Absorption capacity 424%
- Virtually dust free
- Conforms to GLP
- Chip size: 2.8-1.0mm
- Autoclavable

Available in:

50L Polypropylene sack

Product Code:

ECO7-D

Cage Fill Calculator

Cage Size (cm²)

		500 cm ²	1000 cm ²	1575 cm ²	2000 cm ²
Bedding Depth	1 cm	Up to 150 cages	75	48	37
	1.5 cm	100	50	31	25
	2 cm	75	37	24	19



eco-pure™

the natural choice



www.datesand.com

T: 0844 579 0999 (UK local rate)
T: +44 (0)161 274 1080 (International)

