

mini-thermacage⁺

a datesand product

small animal warming cabinet



Supplied by



Allied Scientific
PRODUCTS

www.alliedscientific.com.au

datesand
group

www.datesand.com



Introducing the improved mini-thermacage

The advanced animal warming system with the smaller footprint.



Thermacage Features

Small footprint

320mm W, 400mm L, only 150mm H (300mm H with heat unit).

Integrated

P.I.D. controller, easy adjustment and 3 step operation.

Accurate

precise temperature calibration with test certification.

Easy access

high front, spring loaded, hinged access doors for animal security.

Removable

the dividers provide 4, 3 or 2 compartment configurations.

Easy clean

removable heat unit makes cleaning and relocation easy.

Transparent

clear easy-view end panels and animal friendly translucent dividers.

Reliable

simple, efficient and robust design with 12 months manufacturer's warranty.

Specification

180w heat unit, plug and 3m mains lead. PAT (Portable Appliance Tested). Conforms with Health & Safety and Electricity at Work regulations and ISO standards. 1-100°C temperature scale accurate to +/- 3%. Calibrated/certified with 12 month warranty. CE tested and marked.



This top quality warming cabinet has been designed for use with IVC units where space is at a premium. Its precision temperature control, simple to use 3 step instructions and easy animal access, make it the perfect solution to post-operative care and blood taking procedures.


Now with additional features for better animal welfare and user operation.

Also available, the Thermacage Standard, designed for larger rodent recovery procedures, where space-saving is not as critical. It offers similar specifications, with much larger convalescing compartments.



New Features

1




New diffuser for better airflow and even heat distribution



2

Lockable lid secures 2, 3 or 4 sections



3

New dampers reduce closure noise by up to 15dB